

# OPTIONS FOR A COMPLETE PACKAGE

## VACUUM PUMPS

Totally oil free, mechanical (diaphragm) vacuum pumps with features of advanced design, low noise level, high efficiency and long lifespan.

The LabTech vacuum pumps are specially designed for laboratory operations to meet the highest expectations of precision, reliability and ease of use. They can be used in a multitude of applications, such as vacuum filtration, degassing, impregnation, drying chambers, rotary evaporations, distillation, gel dryers, gas analytical systems and roughing of wide-range turbo molecular pumps.



	VP30	VP18 PLUS	VP50 PLUS
Power (W)	80	30	100
Pumping speed (L/min)	30	18	50
Vacuum pressure (mbar/psi)	250/3.62	20/0.29	150/2.17
Vacuum controller	-	Manual	Manual
Noise (dB)	≤50	≤42	≤58
Dimensions (LxWxH cm)	30x17x24	30x17x24	35x18x24
Weight (Kg)	3	3	6
Power supply voltage	115V 50/60Hz		
Ordering information	LV5010	LV5080	LV5050